


accuslit▶

ophthalmic surgical blades

Accuslit. ophthalmic surgical blades are designed to deliver exceptional precision, sharpness and consistency during ophthalmic and specialty procedures.

Eyegen Medical is proud to offer a broad range of conventional blades for cataract surgery. Keratome/Slit. Sideport. Crescent and MVR blades come in a wide range of sizes to meet the needs of ophthalmic surgeons.



Keratome/Slit Blades



Accuslit Keratome/Slit blades, available in both single & double bevel designed to create accurate width incision for easy/trouble-free penetration. 10 per box.

Single Bevel	Double Bevel	Description
K18011	K18012	1.8 mm, 45° Angled, Bevel Up
K22011	K22012	2.2 mm, 45° Angled, Bevel Up
K25011	K25012	2.5 mm, 45° Angled, Bevel Up
K26011	K26012	2.6 mm, 45° Angled, Bevel Up

Single Bevel	Double Bevel	Description
K27511	K27512	2.75 mm, 45° Angled, Bevel Up
K28011	K28012	2.8 mm, 45° Angled, Bevel Up
K30011	K30012	3.0 mm, 45° Angled, Bevel Up
K32011	K32012	3.2 mm, 45° Angled, Bevel Up

Crescent Blades



Accuslit Crescent blades, angled, offers the flexibility to create optimal tunnel or flap to promote a self-sealing wound. 10 per box.

Item #	Description
C25021	2.5 mm, 55° Angled, Bevel Up
C25023	2.5 mm, 45° Angled, Bevel Down

Sideport Blades



Accuslit Stab / Sideport blades, designed to create precise and true point for easy penetration during initial stab incisions. 10 per box.

Item #	Description
L15002	15°, Straight
L30002	30°, Straight

MVR Blades



Accuslit Vitrectomy/MVR blades, creates precise and fluid tight stab incisions for probe and other posterior segment instruments. 10 per box.

Straight	Angled	Description
M1900	M1901	19G, MVR Blade
M2000	M2001	20G, MVR Blade
M2400	M2401	24G, MVR Blade

Extension Blades (IOL)



Accuslit standard implant blades are designed to enlarge the initial incision for IOL insertion after removal of the crystalline lens. Sides are sharp and tip is blunt. 10 per box.

Item #	Description
K52011	5.2 mm, 45° Angled, Bevel Up

Eye Drape

Eyegen Medical offers EyeProtect drapes with an attached fluid collection bag to remove blood and fluid from the surgical area. Non-woven drapes are lightweight and breathable for optimum patient comfort. They are low linting and moisture resistant to keep patient dry.


Types

- Non-Woven SMS Drape with Fluid Collection Bag
- Non-Woven SMS Drape with Fluid Collection Bag with Pre-Slit

Half Body Drapes

Product Code	Overall size	Adhesive area	Fluid Collection bag size
EM/DB/5040	50 × 40 cm	8 × 6 cm	16 × 12 cm
EM/DB/6060	60 × 60 cm	8 × 6 cm	16 × 12 cm
EM/DB/8080	80 × 80 cm	8 × 6 cm	16 × 12 cm
EM/DB/9090	90 × 90 cm	8 × 6 cm	16 × 12 cm
EM/DB/10070	100 × 70 cm	8 × 6 cm	16 × 12 cm
EM/DB/100100	100 × 100 cm	8 × 6 cm	16 × 12 cm
EM/DB/100120	100 × 120 cm	8 × 6 cm	16 × 12 cm


*Customized sizes are available based on requirement



Full Body Drapes

EyeProtect Full body drapes features wider widths and a two-sided collection bag for maximum fluid capacity. Full body drapes can be utilized in all ophthalmic procedures.

Product Code : EM/FDB/100270
 Overall Size : 100 × 270 cm
 Adhesive Area : 8 × 6 cm
 Collection bag size : 16 × 12 cm




Trolley Drape / Trolley Sheet

EyeProtect Trolley sheets are made of Non-Woven SMS/Poly material. These are convenient to use as they are packed individually in pouch and sterile.

Product Code	Drape size	Material
EM/TB/5075P/S	50 × 75 cm	Non-Woven SMS / Poly
EM/TB/10075P/S	100 × 75 cm	Non-Woven SMS / Poly
EM/TB/100100P/S	100 × 100 cm	Non-Woven SMS / Poly
EM/TB/100120P/S	100 × 120 cm	Non-Woven SMS / Poly
EM/TB/150150P/S	150 × 150 cm	Non-Woven SMS / Poly

*Customized sizes are available based on requirement



Sterile Eye Dressing Pad

Sterile Eye Pads are designed for comfort, absorbency and retention of blood or other discharge when they are applied as postoperative dressing following ocular surgery. Eye Pad is made of White Surgical Cotton with Top Layer highly absorbent spun lace material.

Eye pads are available as:

1. Non - Adherent pad
2. Eye Pad with Adhesive Patch

Packaging: 1 eye pad / sterile pouch; 20 pouches / pack



Eye Shield

Universal Eye Shield

Transparent & coloured, anatomically-shaped eye shield, Features smooth edges and ventilation holes for patient comfort. For use on left or right eye, to protect the eye after surgery. Ready for use for cataract surgery.

Packaging - 10 Shields in a packet



Surgeon Gown

We are engaged in offering an exclusive range of EyeProtect Surgical Gown, designed from the excellent quality Non-Woven SMS material which assures you perfect protection of a surgeon during major surgeries penetration for a short time. Offered surgical gown is a tie-round type, which enables to keep the sterility of surgeons back side. Also, this surgical gown is accessible in various specifications to choose from.

Features:

Tear resistance | Impeccable quality | Long life | Fully Impervious

Product Code	Size
EM/SG/L	Small
EM/SG/XL	Medium
EM/SG/XXL	Large



Patient Gown

EyeProtect Patient Gowns are light weight, breathable and comfortable. It covers body and clothing, and create physical barrier to the transfer of micro organisms and other materials.

Product Code	Size
EM/PG/L	Small
EM/PG/XL	Medium
EM/PG/XXL	Large



Surgeon Cap, Face Mask, Shoe Cover

- i. Surgeon Cap (Elastic & Tie)
- ii. Nurse Cap
- iii. Face Mask (Elastic & Tie)
- iv. Shoe Cover



Eyegen Medical offers an extensive range of eye spears that are designed to meet the needs of clinicians for a variety of ophthalmic procedures. Triangle shaped spearhead, set on a malleable, polypropylene handle. Latex free.

PVA Sponge Spears

High quality PVA Sponge Spears specially designed for lasik and other delicate ophthalmic procedures

Made from exclusive Polyvinyl Acetal (PVA) material. The PVA material is highly absorbent, fast wicking and particulate free. Designed to provide the surgeon with superior fibre free products for unparalleled fluid control.

Packaging - Box containing 250 spears packed in 50 sterile pouches, each pouch containing 5 PVA spears

Product Code: EM/FC001




Cellulose Spears

High quality Cellulose Sponge Spears specially designed for fast and high capacity fluid absorption during ophthalmic procedures

The standard in eye surgical fluid control. Made from highly absorbent, natural, cellulose material. Designed for use in delicate, surgical areas. Ideal replacement for cotton swabs. Designed to maintain rigidity during wicking process. Green handle.

Packaging - Box containing 250 spears packed in 50 sterile pouches, each pouch containing 5 PVA spears

Product Code: EM/FC002



EyeProtect surgery packs are designed and developed in order to maintain and ensure sterility and quality of the highest standards. Our surgery packs are customisable based on the requirement to ease the work flow procedure in the operation theatre

Pre-Op Patient Pack

Content*:

Qty	Description
1	Patient Gown
1	Shoe Cover
1	Patient Cap
1	Patient Towel
1	Patient Wrist Band

Product Code: EM/DK/0001

*Customizable based on requirement

Surgeon Pack

Content*:

Qty	Description
1	Surgeon Gown
1	Shoe Cover
1	Hand Towel
1	Mask
1	Surgeon Cap

Product Code: EM/DK/0002

*Customizable based on requirement

Cataract Pack

Content*:

Qty	Description
1	Trolley drape 150 x 150 cm
1	Eye drape 100 x 120 cm
2	1 ml Syringe
1	5 ml Syringe
1	10 ml Syringe
1	Viscoelastic Cannula
1	Irrigating Cystotome
1	Anterior Chamber Cannula
1	Disposable Simcoe I/A Cannula
1	27G Needle
1	Eye Pad
1	Eye Shield
1	PVA Spears (Pack of 5)
1	Gauze Pads 10 x 10 (Pack of 3)
1	Fluid Administration Set (with Filter)

Product Code: EM/DK/0003

*Customizable based on requirement

IVT Pack

Content*:

Qty	Description
1	Braunstein Caliper
1	Wire Speculum
1	Eye Drape 60 x 60 cm
3	Double ended Cotton buds
1	Eye Pad
1	5 ml Syringe (with needle)
1	1 ml Syringe (with needle)
1	30G IVT needle
2	Sterile gauze
1	Gallipot 60 ml
1	Plastic Tray
1	Wrap Towel 50 x 40 cm

Product Code: EM/DK/0004

*Customizable based on requirement

Lasik Pack

Content*:

Qty	Description
1	Surgery Drape 60 x 40 cm
3	Sterile gauze 5 x 5, 4-ply
5	PVA Spear
1	Eye Shield

Product Code: EM/DK/0005

*Customizable based on requirement

Chalazion Pack

Content*:

Qty	Description
1	Trolley Cover 90 x 90 cm
1	Surgical drape 40 x 40 cm
1	Surgical Tray
1	Syringe 2 ml luer lock
1	Ramplays plastic forceps 14 cm
2	Hand towel 30 x 40
2	Eye pad
2	Gallipots 120 ml
5	Gauzes non-woven 10 x 10 cm
5	Cotton Swabs

Product Code: EM/DK/0006

*Customizable based on requirement