

### Eye Drape

Eyegen Medical offers EyeProtect drapes with an attached fluid collection bag to remove blood and fluid from the surgical area. Non-woven drapes are lightweight and breathable for optimum patient comfort. They are low linting and moisture resistant to keep patient dry.


#### Types

- a. Non-Woven SMS Drape with Fluid Collection Bag
- b. Non-Woven SMS Drape with Fluid Collection Bag with Pre-Slit

### Half Body Drapes

Product Code	Overall size	Adhesive area	Fluid Collection bag size
EM/DB/5040	50 × 40 cm	8 × 6 cm	16 × 12 cm
EM/DB/6060	60 × 60 cm	8 × 6 cm	16 × 12 cm
EM/DB/8080	80 × 80 cm	8 × 6 cm	16 × 12 cm
EM/DB/9090	90 × 90 cm	8 × 6 cm	16 × 12 cm
EM/DB/10070	100 × 70 cm	8 × 6 cm	16 × 12 cm
EM/DB/100100	100 × 100 cm	8 × 6 cm	16 × 12 cm
EM/DB/100120	100 × 120 cm	8 × 6 cm	16 × 12 cm


\*Customized sizes are available based on requirement



### Full Body Drapes

EyeProtect Full body drapes features wider widths and a two-sided collection bag for maximum fluid capacity. Full body drapes can be utilized in all ophthalmic procedures.

Product Code : EM/FDB/100270  
 Overall Size : 100 × 270 cm  
 Adhesive Area : 8 × 6 cm  
 Collection bag size : 16 × 12 cm



### Trolley Drape / Trolley Sheet

EyeProtect Trolley sheets are made of Non-Woven SMS/Poly material. These are convenient to use as they are packed individually in pouch and sterile.

Product Code	Drape size	Material
EM/TB/5075P/S	50 × 75 cm	Non-Woven SMS / Poly
EM/TB/10075P/S	100 × 75 cm	Non-Woven SMS / Poly
EM/TB/100100P/S	100 × 100 cm	Non-Woven SMS / Poly
EM/TB/100120P/S	100 × 120 cm	Non-Woven SMS / Poly
EM/TB/150150P/S	150 × 150 cm	Non-Woven SMS / Poly

\*Customized sizes are available based on requirement

